


CytoFoam[®]

CORE

'FNA — Process Sequence'



A: The FNA foam residue device, CytoFoam Core, consists of a core of absorbent foam housed in a Luer type plastic adapter.

B & C:

A needle is attached to one end and, if a suction FNA technique is to be used, a syringe may be attached to the other end.

D: During the FNA procedure the assembly should be held by the needle hub, or should be connected to a syringe if a suction FNA technique is to be used.

E: The sample is absorbed into the tip of the foam core. Once the FNA sample has been collected the blue Luer type plastic housing should be separated from the needle.

F: The foam core **must** be protected by fitting the plastic guard cap (supplied) over the tip before formalin fixing for at least 6 hours — refer to user instructions for full requirements.

G: After formalin fixation the core is pulled from the adapter, wrapped in processing paper, paraffin processed and sectioned in the usual way.

H: An H&E stained section showing tumour cells within the CytoFoam Core.

Refer to the user instructions for full details



Exmoor Innovations Ltd

Lisieux Way, Taunton, TA1 2LB, UK

T 44 (0) 1823 354307 F 44 (0) 1823 354035

E sales@exmoorinnovations.com

www.exmoorinnovations.com

better diagnosis by design